

Product specification

Cellutec Pro / Pro + · Vacuum tissue processor

Cellutec Pro Item No.: 11000700
 Cellutec Pro + Item No.: 11000710



Specifications

Type	Fast vacuum tissue processor
Principle	Pressure / vacuum
Reagent tanks	12 pcs.
Reagent tanks capacity	5 l
Paraffin tanks	3 pcs.
Paraffin tanks capacity	5 l
Processing capacity	up to 420 standard cassettes
Remote fill & drain ports	2
Wax dumping ports	1
Processing programs	18
Purge reagents	2
Reagents temperature range	20 - 65 °C
Paraffin temperature range	55 - 65 °C
Process start delay time	up to 2 weeks
Reagent partial load capacity	yes
Reagent management system (RMS)	yes
Remote fill & drain (RFD)	yes
External wax dumping (EWD)	yes
USB ports	5

Specifications

Remote alarm ports	1
Network ports (Ethernet)	2
Approvals	CE / IVD / UL
Available software languages	CN, CZ, EN, FR, DE, HU, IT, KZ, PL, PT, RO, RU, ES, TR
Electrical parameters	115 - 230 V, 50/60 Hz
Dimensions (w x d x h), with monitor	680 x 700 x 1500 mm
Weight w/o accessories	226 kg

Standard components

	Pro	Pro +
Sample processing chamber	✓	✓
Sample baskets (3 pcs.)	✓	✓
Reagent tanks 5 l capacity (12 pcs.)	✓	✓
Paraffin tanks 5 l capacity (3 pcs.)	✓	✓
EHE for fast processing		✓
Backlit bottles		✓
NFC (Near Field Communication) technology		✓
LIS integration capability	✓	✓
SPC level sensors (3 pcs.)	✓	✓
Charcoal filter	✓	✓
Color touch screen 15"	✓	✓
Mains cable	✓	✓
Manual	✓	✓

Optional accessories

	Item No.
RFID reader for user identification	11000720
ALIS (connection with LIS system Dedalus Armonia)	11000721
OLIS (connection with other LIS systems)	11000722
UPS (uninterruptable power supply)	11000723
Barcode Reader Gun	11000724
More on request	

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
 55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
 E: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)