

Product specification

Coverttec · Film coverslipper

Coverttec 230 V
Coverttec 115 V

Item No.: 19000000
Item No.: 19000001



Specifications

	Coverttec 230 V	Coverttec 115 V
Nominal power	230 V ~	115 V ~
Nominal frequency	50/60 Hz	60 Hz
Power	120 VA	
Fuse	T 1,6 A	T 2,5 A
Protective class	I	
Overvoltage category	II	
Pollution class	II	
Coverslipping process	up to 1,200 slides per hour	
Coverslipping time	approx. 3 seconds per slide	
Liquids	Xylene	
Manual Xylene dispenser	0.05 – 1 ml	
Coverfilm	max. width 24 mm	
Basket capacity	30 slides	
Filter	Activated carbon filter	
Fan speed	max. 230 m ³ /h	
Film-application length	40 – 60 mm	
Connection	SD Slot for software update	
Dimension (w x d x h)	730 x 310 x 545 mm	
Weight	65.5 kg	

Standard components

Basket (plastic) for 30 standard slides, without hanger, 3 pcs.

Xylene reservoir, 500 ml, 1 pcs.

Xylene waste bottle, 100 ml, 1 pcs.

Film blade – pre-installed, 1 pcs.

Film blade, 1 pack à 5 pcs.

Active carbon filter – pre-installed, 1 pcs.

Coverfilm, 1 roll of 66 m

Mains cable

Manual

Optional accessories

Item No.

Basket (plastic) for 30 standard slides, without hanger, compatible with Cromatec, 1 pcs.

38501051

Basket (plastic) for 30 standard slides, with removable hanger, compatible with Cromatec, 1 pcs.

41000420

Xylene waste bottle, 100 ml

19000102

Xylene bottle, 500 ml

19000103

Consumables

Item No.

Film blade, 1 pack à 5 pcs.

19000100

Activated carbon filter, 1 pcs.

19000101

Coverfilm, 1 roll of 66 m

19000104

Xylene, 5 l

20011070.05000

Xylene, 25 l

20011070.25000

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
E: mail@slee.de



www.slee.de



YouTube



LinkedIn



Instagram