

**TANNER**  
scientific®

# TN50 & TN50UV CRYOSTATS



**IN-STOCK & READY-TO-SHIP!**

## Product Description:

The German-made **Tanner Scientific® TN50 Cryostat** features a smooth microtome for easy, intuitive operation. For extra convenience the system includes a spacious stainless steel chamber, 24 easily accessible freezing positions and a well illuminated cooling chamber that can be effectively cooled down to  $-35^{\circ}\text{C}$ . For additional user safety the system can optionally be equipped with an efficient UV-C disinfection system.



**2 YEAR WARRANTY\***

MADE IN GERMANY

## Item # TNR TN50

**Dimensions (WxDxH):** 26 $\frac{3}{4}$ " x 23 $\frac{5}{8}$ " x 43 $\frac{3}{4}$ "

**Weight (w/o accessories):** 280 lbs

**Electrical Parameters:** 115V; 50/60 Hz

- **Exceptional German engineering**
- Ideal for Mohs surgery and traditional pathology applications
- 2 Year Warranty\*
- Optional UV-C disinfection solution
- Removable Flat Mounting Surface Accessory
- Smooth rotary microtome
- Spacious cooling chamber for efficient and safe workflow
- Glass anti-roll plate
- 24 easily accessible freezing positions
- Heated viewing window
- Easy-to-reach electronic control panel
- Automatic defrost
- Removable condensate bottle with level sensor

**MERCEDES  
SCIENTIFIC**



MercedesScientific.com



800.331.2716



800.359.8807

All products are subject to change, availability or discontinuation at any time without prior notice. All trademarks and brands are the property of their respective owners.