

EZ-Writer

MANUAL PRINTER FOR CASSETTES

MADE IN GERMANY 

item # VOG CP6001

Product Description:

Compact cassette printer for on-demand printing, the unit offers flexibility, consistency and reliability. Can be used with computer, bar code scanner, LAN-access and USB Windows-Driver. Customer has the possibility to use the printer as an online job printer or as a stand-alone unit.



Cassette Printer

2 lines, default font, without bar code:

2 x 16 = 32 characters

2 lines, default font, with bar code:

2 x 11 = 22 characters

1 line, uppercase, without bar code:

1 x 13 = 13 characters

Similarity

- On-demand printing with or without Matrix Code
- Accepts most standard size cassettes
- Equipped with number key PAD for quick entry
- Special characters available
- Recall saved screens from memory
- Stored computer information can be downloaded
- Upper and lower case letters
- Different printing mode
- Capable of printing special characters
- Stores up to 100 different programmable formats
- Bar code scanner adaptable (Y-gate)
- Prints cassette in approximately 2-3 seconds
- LAN-access makes it easy to transfer external data from a computer/laptop to the printer
- USB Windows Driver available
- Ribbon with UV-sensitive substance and RFID-chip
- Size: 25 x 18 x 19 cm



PHONE: 800.331.2716 | FAX: 800.359.8807 | WEB: MERCEDESSCIENTIFIC.COM