

Product specification

Aquatec · Transfer system for microtomes of CUT series

Item No.: 10090056



Specifications

Heated waterbath	up to +50 °C
Volume of heated waterbath	750 ml
Volume of pump reservoir	570 ml
Required liquid	demineralized or distilled water
Flow rate	0 - 500 ml / min.
Operating temperature range	+10 to +35 °C
Operating humidity	max. rel. 80 % non-condensing
Storage temperature range	+5 to +55 °C
Storage humidity	max. rel. 80 % non-condensing
Pollution degree	2
Overvoltage category	II
Power requirements	100 - 240 V / 50/60 Hz, 0.7 A + / -10 %
Weight	
basic device w/o transfer bridge, blade holder, incl. lid	5.0 kg
disposable blade holder AQ	2.4 kg
Dimensions (W x D x H) w/o transfer bridge, blade holder	
basic device mounted on microtome	300 x 235 x 128 mm
basic device unmounted on microtome	300 x 280 x 128 mm

Standard components

Basic device incl. removable tray

Disposable blade holder AQ

Lightning

Waste tray

Holder for temperature sensor

Temperature sensor

2 x transfer inserts (1 x installed)

2 x suction filter (1 x installed)

Syringe

Drain screen

Manual

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
M: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)